



HbA1c Fast Test Kit (Immunofluorescence Assay)

User Manual

Getein1100: Cat.# IF1017
Getein1600: Cat.# IF2017

INTENDED USE

HbA1c fast test kit is intended for the quantitative measurement of HbA1c in whole blood. This test is used as an aid for monitoring glycemic control in diabetics. In addition, it can identify people at risk of developing the disease and ongoing monitoring.

SUMMARY

Hemoglobin is the protein molecule in red blood cells with the main function transport oxygen and carbon dioxide in blood. HbA1c belongs to the glycosylated hemoglobin, a fraction formed by the attachment of various sugars to the Hb molecule and is proportional to average blood glucose concentration over the previous four weeks to three months. One advantage of using HbA1c for diagnosis is that the test does not require a fasting blood sample. Although HbA1c testing is mainly used for monitoring blood sugar control in patients with diabetes, the World Health Organization (WHO) now recommends that HbA1c can be used as a diagnostic test for diabetes, provided that stringent quality assurance tests are in place and assays are standardised to criteria aligned to the international reference values.

PRINCIPLE

The test uses an anti-human Hb monoclonal antibody conjugated with fluorescence latex and an anti-human HbA1c monoclonal antibody coated on the test line. After the sample has been applied to the test strip, the fluorescence

latex-labelled anti-human Hb monoclonal antibody binds with the HbA1c and Hb in sample proportionally and forms marked antigen-antibody complex. The complex moves to the detection zone by capillary action. Then marked antigen-antibody complex is captured on the test line by the anti-human HbA1c monoclonal antibody. The fluorescence intensity of the test line increases in proportion to the amount of HbA1c in sample. Then insert test card into Getein1100 Immunofluorescence Quantitative Analyzer/Getein1600 Immunofluorescence Quantitative Analyzer (hereinafter referred to as Getein1100 and Getein1600), the concentration of HbA1c in sample will be measured and displayed on the screen. The value will be stored in Getein1100/Getein1600 and available for downloading. The result can be easily transmitted to the laboratory or hospital information system.

CONTENTS

- A kit for Getein1100 contains:
 - Getein HbA1c test card in a sealed pouch with desiccant 25
 - Disposable pipet 25
 - Sample diluent 25
 - SD card 1
 - User manual 1
- A kit for Getein1600 contains:
 - Sealed cartridge with 24/48 Getein HbA1c test cards 2
 - User manual 1
 - Package specifications:
 - 2×24 tests/kit, 2×48 tests/kit
 - Materials required for Getein1600:
 - Sample diluent 1
 - Box with pipette tips 1
 - Coated wells 1
- Sample diluent composition:
 - Phosphate buffered saline, proteins, detergent, preservative, stabilizer.
- A test card consists of:

A plastic shell and a reagent strip which is composed of a sample pad, nitrocellulose membrane (one end of the membrane is coated with a fluorescence latex-labelled anti-human Hb monoclonal antibody, the test line is coated with an anti-human HbA1c monoclonal antibody, and the control line is

coated with rabbit anti-mouse IgG antibody), absorbent paper and liner.

Note: Do not mix or interchange different batches of kits.

APPLICABLE DEVICE

Getein1100 Immunofluorescence Quantitative Analyzer
Getein1600 Immunofluorescence Quantitative Analyzer

STORAGE AND STABILITY

Store the test card at 4–30°C with a valid period of 24 months. Use the test card for Getein1100 within 1 hour once the foil pouch is opened.

Use the test card for Getein1600 within 24 hours once opened. Store the sample diluent at 0–30°C with a valid period of 24 months.

Store the sample diluent at 2–8°C for better results.

PRECAUTIONS

- For *in vitro* diagnostic use only.
- For professional use only.
- Do not use the kit beyond the expiration date.
- Do not use the test card if the foil pouch or the cartridge is damaged.
- Do not open pouches or the cartridge until ready to perform the test.
- Do not reuse the test card.
- Do not reuse the pipet.
- Handle all specimens as potentially infectious. Proper handling and disposal methods should be followed in accordance with local regulations.
- Carefully read and follow user manual to ensure proper test performance.

SPECIMEN COLLECTION AND PREPARATION

- This test can be used for *whole blood samples*. *Heparin, sodium citrate and EDTA* can be used as the anticoagulant under aseptic conditions.

- The test is for human blood, other specimens or bodily fluids may not get accurate results.
- The test should be performed within 4 hours after whole blood collection.
- Samples could be kept for 7 days at 2~8°C and avoid cryopreservation.
- Samples must be recovered to room temperature before testing.
- SAMPLE VOLUME (for Getein1100): 10 µl.**

TEST PROCEDURE

- Collect specimens according to user manual.
- Test card, sample should be brought to room temperature before testing.
- For Getein1100:**
 - Confirm SD card lot No. in accordance with test kit lot No.. Perform "QC (SD)" calibration when necessary (Details refer to 8.5.2 of Getein1100 User Manual).
 - On the main interface of Getein1100, press "ENT" button to enter testing interface.
 - Remove the test card from the sealed pouch immediately before use. Label the test card with patient or control identification.
 - Put the test card on a clean table, horizontally placed.
 - Using sample transfer pipette, deliver 10 µl of sample into one tube of sample diluent, mix gently and thoroughly. Then drop 100 µl of sample mixture (or 3~4 drops of sample mixture when using disposable pipet) into the sample port on the test card.
 - Reaction time: 5 minutes.** Insert the test card into Getein1100 and press "ENT" button after reaction time is elapsed. The result will be shown on the screen and printed automatically.
- For Getein1600:**
 - Each cartridge for Getein1600 contains a specific RFID card which can calibrate automatically.
 - Place samples in the designed area of the sample holder, insert the holder and select the right test item, Getein1600 will do the testing and print the result automatically.

Notes:

- It is required to perform "SD Card Calib" calibration when using a new batch of kits.

- It is suggested to calibrate once for one batch of kits for Getein1100.
- Make sure the test card and the sample insertion is correct and complete.

TEST RESULTS

Getein1100/Getein1600 can scan the test card automatically and display the result on the screen. For additional information, please refer to the user manual of Getein1100/Getein1600.

EXPECTED RANGE OF VALUE

HbA1c concentration is determined using samples obtained from 345 apparently healthy individuals. The normal value for HbA1c is 3.8%-5.8%. It is recommended that each laboratory establish its own expected values for the population it serves.

PERFORMANCE CHARACTERISTICS

Measuring Range 2%-14%
 Lower Detection Limit ≤2%
 Within-Run Precision (n=10) ≤10%
 Between-Run Precision ≤15%
 Accuracy: verify with comparison experiments, the correlation coefficient $r \geq 0.990$, the relative error ≤20%.

LIMITATIONS

- The result of the test should be evaluated in the context of all the clinical and laboratory data available. In those instances where the laboratory results do not agree with the clinical evaluation, additional tests should be performed accordingly.
- Some substances in blood as listed below may interfere with the test and cause erroneous results. The maximum allowance concentration of them is as follows:

Interferent	Concentration (Max)
Triglyceride	25 g/L
Bilirubin	0.1 g/L

REFERENCES

- Cagliero E, Levina E V, Nathan D M. Immediate feedback of HbA1c levels improves glycemic control in type 1 and insulin-treated type 2 diabetic patients[J]. *Diabetes care*, 1999, 22(11): 1785-1789.
- Özdamar Ö, Gün i, Keskin U, et al. The role of maternal serumbeta-HbA1c and PAPP-A levels at gestational weeks 10 to 14 in the prediction of pre-eclampsia[J]. 2014.

DESCRIPTION OF SYMBOLS USED

The following graphical symbols used in or found on HbA1c Fast Test Kit (Immunofluorescence Assay) are the most common ones appearing on medical devices and their packaging. They are explained in more details in the European Standard EN 980:2008 and International Standard ISO 15223-1:2007.

Key to symbols used			
	Manufacturer		Expiration date
	Do not reuse		Date of manufacture
	Consult instructions for use		Batch code
	Temperature limitation		In vitro diagnostic medical device
	Sufficient for		Authorized representative in the European Community
	CE mark		Do not use if package is damaged

Thank you for purchasing HbA1c Fast Test Kit (Immunofluorescence Assay). Please read this user manual carefully before operating to ensure proper use.

Version: WIF22-S-01

Getein Biotech, Inc.
 Add: No.9 Bofu Road, Luhe District, Nanjing, 211505, China
 Tel: +86-25-68568508
 Fax: +86-25-68568500
 E-mail: tech@getein.com.cn
 overseas@getein.com.cn
 Website: www.bio-GP.com.cn